

Accuracy of the Anti-SARS-VoV-2 Rapid Test Kit

All testing was compared to RT-PCR tests performed from naso- or oropharyngeal swabs specimens.

Table 8: SARS-CoV-2 **IgM** Positive results of three disease periods from symptom onset

Disease period (days)	Samples No. of PCR positive	No. of positive results by Autobio	Clinical Sensitivity	95%CI
≤7	51	19	37.25%	23.99%~50.52%
8-14	52	38	73.08%	61.02%~85.13%
≥15	302	289	95.70%	93.41%~97.98%

Table 9: SARS-CoV-2 **IgG** Positive results of three disease periods from symptom onset

Disease period (days)	Samples No. of PCR positive	No. of positive results by Autobio	Clinical Sensitivity	95%CI
≤7	51	16	31.37%	18.64%~44.11%
8-14	52	34	65.38%	52.45%~78.32%
≥15	302	299	99.01%	97.89%~100.13%

Table 6: Results of **IgM** Clinical Sensitivity and Specificity – for all time periods from symptom onset

Evaluation reagent		PCR		Total
		Positive	Negative	
Autobio	Positive	346	1	347
	Negative	59	311	370
Total		405	312	717
IgM Clinical Sensitivity		85.43% (95%CI: 82.00%-88.87%)		
IgM Clinical Specificity		99.68% (95%CI: 99.05%-100.31%)		

Table 7: Results of **IgG** Clinical Sensitivity and Specificity – for all time periods from symptom onset

Evaluation reagent		PCR		Total
		Positive	Negative	
Autobio	Positive	349	2	351
	Negative	56	310	366
Total		405	312	717
IgG Clinical Sensitivity		86.17% (95%CI: 82.81%-89.53%)		
IgG Clinical Specificity		99.36% (95%CI: 98.47%-100.24%)		